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1023

Revision History

Revised Chapter	Content Revised	Date of Revision	Revised By	Number of Revision	Remarks

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MENU

Using Range of This Instruction	1
Information About This Autoclave	1
Explanation Of Symbols	2
Safety Warning	3
1. Brief Introduction	4
2. Using Range	4
3. Parameters and Sterilization Programs	5
4. Operation Panel	6
5. Installation	13
6. Operation	14
7. Abnormal Situations	16
8. Maintenance	17
9. Transportation And Storage	20
10. Warranty	21
11. Accessories	22
Appendixes	23
Appendix 1: Structure Diagram	23
Appendix 2: Circuit Diagram	24
Appendix 3: EMC	25

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	Symbol for "Caution!"
	Symbol for "PROTECTIVE CONDUCTOR TERMINAL"
	Symbol for "HOT SURFACE, KEEP AWAY"
	Symbol for "ENVIRONMENTAL PROTECTION – Electrical wastes should not be disposed together with household wastes. If possible, please contact local distributor or government to recycle it."
	Symbol for "MANUFACTURER"
	Symbol for "COMPILES WITH MDD93/42/EEC REQUIREMENTS"
	Symbol for "DATE OF PRODUCTION"
	Symbol for "SERIAL NUMBER"
	Symbol for "EUROPEAN REPRESENTATIVE"
	Symbol for "THIS SIDE UP"
	Symbol for "KEEP AWAY FROM RAIN"
	Symbol for "DO NOT ROLL"
	Symbol for "STACKING LIMITED 3"
	Symbol for "Room Temperature: 5°C~40°C"
	Symbol for "Relative humidity: ≤80% "

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Electric shock, fire or equipment damage might happen if you ignore below warnings.

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Type of this autoclave	
Type	Description of intended use
B	The sterilization of all wrapped or non-wrapped, solid, hollow load products type A and porous products as represented by the test loads in standard: EN13060:2004+A2:2010

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
The data have "*" can be adjusted

B&D Test: The display shows 4min, but normally it takes 3.5min.

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Please ensure the table for installation is strong enough for an autoclave.

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The autoclave might not run if the user does not adjust its atmospheric pressure.

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Please ensure the autoclave is well grounded.

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Only distilled water can be used for the autoclave. Otherwise, there will be problem with the steam generator and solenoid valves. The user shall be responsible for the problem if he/she uses other kinds of water.

Before filling water, please turn on the autoclave and ensure the autoclave is horizontal. After filling water, if possible, please drain the waste water.

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If the temperature and pressure inside the chamber are high, the door will be difficult to be closed. We suggest to close the door quickly, or, you can wait 30 seconds, then close the door.

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The user shall fill in water promptly if there is a low-water alert. Otherwise it will show "E08 or E9" error alarm.

You shall use a tray-holder tool when taking out loads to prevent scald. Don't open the door until the pressure is within "-05~05".

To ensure the effectiveness of sterilization, during daily use, we suggest putting test paper or pouches with indicators together with the loads into the sterilization chamber.

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The waste water might be very hot, please be careful when you are draining them.

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We suggest running program once more to see if the error happens again.

No.	Code	Sound	Conditions	Reason	Resolution
1	E31	"Du"long beep	Temperature in chamber >150℃;	The sensor of steam generator does not work.	Check temperature sensor in chamber
2	E32	"Du"long beep	Temperature of Chamber heater>280℃;	The sensor of steam generator does not work.	Check temperature sensor of the chamber heater
3	E51	"Du"long beep	Temperature in chamber ≤0℃;	The sensor of steam generator does not work.The room temperature is too low.	Check temperature sensor in chamber Check room temperature.
4	E52	"Du"long beep	Temperature of chamber heater≤0℃;	The sensor of steam generator does not work.The room temperature is too low.	Check temperature sensor in chamber Check room temperature.
5	E63	"Du"long beep	1、 steam generator temperature≤0℃; 2、 steam generator temperature>230℃;	steam generator temperature>230℃; Steam temperature control instability, over 230℃, steam generator temperature sensor damaged.	Check steam generator temperature sensor, control board, steam generator
6	E2	"Du"long beep	The sterilization pressure is 40KPA more than designed pressure. (134℃:210kpa /121℃process:110kpa)	vacuum pump damaged.	Check vacuum pump. Do a vacuum test
7	E61	"Du"long beep	134℃process: inner temperature >140℃ or 121℃ process: inner temperature >127℃; temperature control instability.	temperature sensor inside chamber damaged.	Check temperature sensor in chamber.
8	E62	"Du"long beep	Temperature of Chamber heater>155℃; temperature control instability, control board damaged.	temperature sensor is damaged. Circuit board is damaged.	Check temperature sensor, chamber heater, control board, chamber heater
9	E41	"Du"long beep	In preheat period, after 8mins warm-up, temperature Chamber heater<100℃; chamber heater damaged.	Chamber heater is damaged	check the chamber heater.
10	E42	"Du"long beep	In preheat period, after 8mins steam generator temperature <110℃; heating rod damaged.	Chamber heater is damaged	Check heating rod
11	E5	"Du"long beep	When the period of "exhaust", after working 10mins, the pressure in chamber still over 0.5bar; air relief instability.	Outlet solenoid valve is blocked.	Check water drain valve
12	E6	"Du"long beep	The door is opened during working; the door detector switch damaged.	Door sensor is damaged	Check door sensor.
13	E7	"Du"long beep	The local air pressure is < 70KPa;	local air pressure is too low. Wrong air pressure.	Can not use in these area. Adjust atmospheric pressure:see 5.2
14	E8	"Du"long beep	In pre-vacuum period, every 5mins temperature raise < 3℃.	Air leakage. No water inside tank.	Check water pump, heating rod, control board.Check water level inside tank.
15	E9	"Du"long beep	In sterilization period, the sterilization pressure is 0.3 bar lower than the designed pressure.	No water inside tank. Heater inside steam generator is damaged. Water pump damaged.	Check heaters, water pump, water level.
16	E10	"Du"long beep	The electronic locker is in wrong condition	Electronic locked damaged. Circuit board damaged.	check electronic locker and circuit board
17	E11	"Du"long beep	The electronic locker is in wrong condition	Electronic locked damaged. Circuit board damaged.	check electronic locker and circuit board
18	E12	"Du"long beep	The vacuum not reach -70Kpa 2 times during the program which have at least 3 times vacuum.	Vacuum pump damaged.Fan damaged.	Check vacuum pump, fan
20	E99	"Du"long beep	The communication between CPU is wrong.	Connect disabled. CPU loosed.	Check control board data line, and CPU

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
Item	Checked by	Cycle	Maintenance
Door	Professional engineer	2 years	See 8.1
Seal ring	User	1 year	See 8.5
Printer papers	User	When then printer has no paper	See 8.2
Fuse	Professional engineer	When the fuse is damaged	See 8.6
Safety valve	Professional engineer	1 year	See 8.9

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
After the door adjustment, it is necessary to do a vacuum test. If there is a leakage, the user shall adjust it again.

If the door is too loose, the sealing ring might be exploded out with a big "bang". So, please be careful to use this function!

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Only one side of the paper can be printed on.

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Check the new fuse if it is suitable for this autoclave before replacement.

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(pic8)

- 1, Removing the part 1 in pic8, then removing the pipe which connects the safety valve.
- 2, Removing the screw (part2 in pic8);
- 3, Replace the new safety valve.

The new safety valve shall be same model. If you need any help, please contact your local distributor or us.

Only maintain the autoclave when it cools down by well trained person.

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11. Accessories

No.	Accessories	Q'ty (PC)
1	Drainage pipe(LF-9-10-1)	1
2	Tray (see table 11.1)	3
3	Cable with plug(LF-9-12-2)	1
4	Tray shelf (see table 11.2)	1
5	Tray-hand-holder(LF-9-12-4)	1
6	Fuse(Φ 6X30 230V20A)	2
7	Manual	1


Table 11.1: parts for different models



Part/Model	LF-12L-E	LF-18L-E	LF-22L-E
Tray	LF-2-6-1	LF-1-6-1	LF-3-6-1
Tray shelf	LF-2-6-2	LF-1-6-2	LF-3-6-2

Appendixes

Appendix 1: Structure diagram

Appendix 2: Circuit diagram

Vacuum pump power: AC: 230V 50Hz 70W
 Water pump power: AC: 230V 50Hz 47W
 Heating circle: AC: 230V 50Hz 1500W
 Heating rod: AC: 230V 50Hz 750W
 Electromagnetic valve: DC: 24V 5W

Appendix 3: EMC

Electromagnetic emissions

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Steam sterilizer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Steam sterilizer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Electromagnetic immunity

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	<input type="checkbox"/> 6 kV contact <input type="checkbox"/> 8 kV air	<input type="checkbox"/> 6 kV contact <input type="checkbox"/> 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	<input type="checkbox"/> 2 kV for power supply lines <input type="checkbox"/> 1 kV for input/output lines	<input type="checkbox"/> 2 kV for power supply lines <input type="checkbox"/> 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	<input type="checkbox"/> 1 kV line(s) to line(s) <input type="checkbox"/> 2 kV line(s) to earth	<input type="checkbox"/> 1 kV line(s) to line(s) <input type="checkbox"/> 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Steam sterilizer requires continued operation during power mains interruptions, it is recommended that the Steam sterilizer be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE *UT* is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity

The Steam sterilizer is intended for use in the electromagnetic environment specified below. The customer or the user of the Steam sterilizer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Steam sterilizer including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1,2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,3 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Steam sterilizer is used exceeds the applicable RF compliance level above, the Steam sterilizer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Steam sterilizer.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Steam sterilizer

The Steam sterilizer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Steam sterilizer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Steam sterilizer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Distributed by:

ZKDENTAL Italia

via De Gasperi 85 – 20017 RHO (MI) Italy

MDS MEDICAL LTD

93 headlands, Kettering, Northamptonshire, NN15 5EG, UK